



STEINDORFF®

NYMC62B00

Polarizing Microscope




NEW YORK MICROSCOPE COMPANY INC.
AKA MEL SOBEL MICROSCOPES

100A Lauman Ln., Hicksville, N.Y. 11801 Toll Free: (877) 877-7274 • Fax: (516) 801-2046
Web Site: www.microscopeinternational.com • www.nyscopes.com • E-mail: Info@nyscopes.com



STEINDORFF®



List of Parts

- | | |
|----------------------------------|---|
| 1. Eyepiece | 11. Brightness Control |
| 2. Binocular Head | 12. Power Switch |
| 3. Bertrand Lens | 13. Power Indicator |
| 4. Compensator Socket | 14. Light Collector |
| 5. Clamp Screw | 15. Graduated Marks of Polarizer |
| 6. Stage Height Adjustment Screw | 16. Iris Diaphragm Lever |
| 7. Coarse Focusing Knobs | 17. Stage Centering Screws (Detachable) |
| 8. Fine Focusing Knobs | 18. Nosepiece |
| 9. Condenser Centering Screws | 19. Graduated Marks on Analyzer |
| 10. Condenser Focusing Knob | 20. Analyzer Lever |

STEINDORFF®

Description

The Polarizing Microscope is one of the most useful professional experiment instrument for use in the fields of geology, mineral and metallurgical department as well as associated educational institutions of higher studies.

This Polarizing Microscope can be used for single polarizing examination, orthogonal polarizing examination and conical polarizing examination as well as for photomicrography.

This Microscope is equipped with various accessories such as a gypsum test plate (wavelength), a mica test plate (quarter wavelength).

Specifications

1. Observation Tube - Binocular Head.
2. Objectives: stress less achromatic objectives.

Magnifying Power Objectives	Numerical Aperture N.A.	Thickness of Cover Glass
4x	0.10	-
10x	0.25	-
20x	0.40	
40x	0.65	0.17
60x(optional)	0.85	0.17

3. Eyepieces

Type	Magnifying Power
Standard Eyepiece (Pair)	10x
Eyepiece with Crosshairs	5x
Eyepiece with Graticule	10x

4. Final magnifying powers of objectives and eyepieces combination.

Final magnifying Powers:

Eyepiece	Objectives			
	4x	10x	40x	60x
10x	4x	10x	40x	60x
Magnification	40x	100x	400x	600x

5. Quadruple centralized revolving nosepiece for objectives.
6. Rotating Stage: 130mm in outside dia., calibrated in 360°, Vernier scale 0.1°.
7. Focusing Mechanism: Coaxial co-guide way coarse and fine focusing adjustment.

Focusing Range: 25mm, Division 0.002mm.

STEINDORFF®

8. Illumination: Transmitted 6V-20W light source, brightness adjustable.
9. Line Power: 220V/50-60Hz.
10. Binocular Viewing Head
11. Compensator: Quarter Wedge 1-4 Order, Full Wave 11 Gypsum & Quarter Wave 1/41 Mica.
12. Condenser: Abbe two-lens condenser, N.A. = 1.25.

Unpacking

Identify the version of the Microscope purchased, unpack and bring the Microscope out of the case. Check up the accessories in the case.


Operation

Insert the plug of the power of this microscope into a suitable grounded 220V AC line power socket. Turn on the power by the switch & the indication light will be lit. Switch ON the Microscope Illuminator.

- 1) Place a specimen onto the stage. Make an observation, pay attention to pull the lever for pulling the analyzer out from the light path. Open the Iris Diaphragm. Make sure that the field Iris Diaphragm coincides with the field of view.

- 2) Secure the specimen properly with the clips on the Rotary Stage.

- 3) Adjust the center of the stage using the stage centering screw. Focus the specimen sharply. Find a marked feature in the field of view and make it to situate at the center of the ocular's cross hairs. Rotate the stage, if the optical axis of the objectives not coincides with the center of the rotation of the stage, then the target selected will rotate about certain center S (that is the rotating center of the stage), the trajectory of which is a circle. Turn the target point to the point O₁ to move toward the point S and coincide with it. Then, turn the stage again, to view whether the two points coincide with each other, if any deviation exists, repeat this procedure again.



STEINDORFF®

- 4) When observing the microstructure of a specimen, usually use a low power objective to look for the object, then move the object towards the center of the field of view and exchange a high power objective for observation. Care should be taken to avoid the objective being knocked against the specimen and damaged. When focusing, lift the stage allow it close the objective lens as near as possible. Lower the object surface while viewing. When an image is observed, use the fine adjustment knob for adjustment until a sharp image is obtained.
- 5) Examining the specimen in orthogonal polarized light. As mentioned in section (6) when a sharp image is obtained, because the polarizer is still inserted in the optical path, actually it is in a single polarized light. Then insert an analyzer. Pay attention to the graduated marks of the polarizer and analyzer. When both marks are situated at "0" position, it means that they are in orthogonal. Now, the polarized orientation of the analyzer is in north - south (or up and down).
- 6) Examining the specimen in the conical polarized light. Usually, conical polarized light is used for high power examination. Insert Bertrand lens into the light path when using orthogonal polarizing light. The swinging upper lenses of the condenser are arranged so that they may be moved into the optical axis of the microscope for examining the conical characteristic of the specimen.
- 7) When an indistinct image of the filament of the bulb reveal in the field of view, add a ground glass properly. When the light is yellowish, a filter can be added. (Place a filter above the collector in the base of the microscope in transmitted light).
- 8) Adjusting the light housing.
When the bulbs is replaced, the light housing should be realigned. Insert the light housing into the adapter of the vertical illuminator. Switch on the power, moves the light housing back & forth to eliminate the image of the filament in the field of view. After operating as mentioned in Section (6) then observe through eyepiece tube with the eyepiece removed, until a bright spot at the back focal plane of the objective or the image of the aperture diaphragm can be seen.
Adjust screws at the rear and side surface for centering the filament of the bulb.

USE OF BERTRAND LENS IN POLARIZING MICROSCOPE

The Bertrand Lens is generally used to observe interference figures in high power objective. Please check the following points while using Bertrand Lens.

STEINDORFF®

- 1) Analyzer & Polarizer shall be crossed in the light path. The images appear in the objective rear focal plane when an optically anisotropic specimen is viewed between crossed polarizer using a high numerical aperture objective/condenser combination.
- 2) The iris diaphragm should be closed to 90%.
- 3) The Bertrand Lens shall be in the center to the plane. Please use centering screws provided with the Bertrand Lens.

MAINTENANCE

1. Do not dismantle any part of the instrument
If there is any dirt on the surfaces of the lens, wipe it off gently with some absorbent cotton moistened with xylene. Dust on the lens surfaces can be blown off by using an air blower.
2. Oil the mechanical moving parts of the instrument with some watch or clock lubricant.
3. Avoid contacting high temperature.
4. Objectives and eyepieces should replace back to their containers and place dust caps and dust plates to all of the light ports after using.
5. The instrument should be kept in a cool, shady and dry place, free from dust, acid alkalis fume or vapors.
6. Check and maintain the microscope at regular intervals.

NEW YORK MICROSCOPE COMPANY INC.
AKA **MEL SOBEL MICROSCOPES**

100A Lauman Ln., Hicksville, N.Y. 11801 Toll Free: (877) 877-7274 • Fax: (516) 801-2046
Web Site: www.microscopeinternational.com • www.nyscopes.com • E-mail: Info@nyscopes.com

